

Hwan Myung Kim

List of Publications by Year in descending order

Source: <https://exaly.com/author-pdf/4365037/publications.pdf>

Version: 2024-02-01

116
papers

8,965
citations

36303

51
h-index

42399

92
g-index

118
all docs

118
docs citations

118
times ranked

7614
citing authors

#	ARTICLE	IF	CITATIONS
1	Observing hepatic steatosis with a commercially viable two-photon fluorogenic probe. <i>Materials Chemistry Frontiers</i> , 2022, 6, 553-560.	5.9	19
2	A coumarin-based reversible two-photon fluorescence probe for imaging glutathione near <i>N</i> -methyl-D-aspartate (NMDA) receptors. <i>Chemical Communications</i> , 2022, 58, 3633-3636.	4.1	11
3	Highly Sensitive Two-Photon Lipid Droplet Tracker for <i>In Vivo</i> Screening of Drug Induced Liver Injury. <i>ACS Sensors</i> , 2022, 7, 1027-1035.	7.8	19
4	An azo dye for photodynamic therapy that is activated selectively by two-photon excitation. <i>Chemical Science</i> , 2021, 12, 427-434.	7.4	33
5	Recent progress in the two-photon fluorescent probes for metal ions. <i>Coordination Chemistry Reviews</i> , 2021, 427, 213574.	18.8	85
6	Development of two-photon fluorescence probe for detecting cyclooxygenase-2 level in human colorectal cancer tissue. <i>Sensors and Actuators B: Chemical</i> , 2021, 330, 129329.	7.8	7
7	Analyzing Nonmelanoma Skin Cancer Using Enzyme-Activatable Two-Photon Probes. <i>Bulletin of the Korean Chemical Society</i> , 2021, 42, 103-106.	1.9	5
8	Two-photon ESIPT-based fluorescent probe using 4-hydroxyisoindoline-1,3-dione for the detection of peroxynitrite. <i>Chemical Communications</i> , 2021, 57, 11084-11087.	4.1	37
9	Two-Photon Fluorescent Probes for Detecting Enzyme Activities in Live Tissues. <i>ACS Applied Bio Materials</i> , 2021, 4, 2957-2973.	4.6	17
10	Hypochlorite-Activated Fluorescence Emission and Antibacterial Activities of Imidazole Derivatives for Biological Applications. <i>Frontiers in Chemistry</i> , 2021, 9, 713078.	3.6	6
11	A Diagnostic Method for Gastric Cancer Using Two-Photon Microscopy With Enzyme-Selective Fluorescent Probes: A Pilot Study. <i>Frontiers in Oncology</i> , 2021, 11, 634219.	2.8	7
12	Fluorescence Probe for Imaging <i>N</i> -Methyl-D-aspartate Receptors and Monitoring GSH Selectively Using Two-Photon Microscopy. <i>Analytical Chemistry</i> , 2021, 93, 11612-11616.	6.5	26
13	A red-emissive two-photon fluorescent probe for mitochondrial sodium ions in live tissue. <i>Chemical Communications</i> , 2021, 57, 8929-8932.	4.1	5
14	Near-Infrared Ratiometric Two-Photon Probe for pH Measurement in Human Stomach Cancer Tissue. <i>ACS Applied Bio Materials</i> , 2021, 4, 2135-2141.	4.6	14
15	Highly selective two-photon fluorescent off-on probes for imaging tyrosinase activity in living cells and tissues. <i>Chemical Communications</i> , 2021, 57, 6911-6914.	4.1	12
16	Azulene-based fluorescent chemosensor for adenosine diphosphate. <i>Chemical Communications</i> , 2021, 57, 10608-10611.	4.1	10
17	Highly Stable Red-Emissive Ratiometric Probe for Monitoring β -Galactosidase Activity Using Fluorescence Microscopy and Flow Cytometry. <i>Analytical Chemistry</i> , 2021, 93, 14778-14783.	6.5	19
18	Cancer-Targeted Azo Dye for Two-Photon Photodynamic Therapy in Human Colon Tissue. <i>Analytical Chemistry</i> , 2021, 93, 16821-16827.	6.5	7

#	ARTICLE	IF	CITATIONS
19	Combining hydrophilic and hydrophobic environment sensitive dyes to detect a wide range of cellular polarity. <i>Chemical Science</i> , 2020, 11, 596-601.	7.4	48
20	Asymmetric cyanine as a far-red fluorescence probe for mitochondrial viscosity. <i>Dyes and Pigments</i> , 2020, 174, 108080.	3.7	39
21	Screening of Drug-Induced Steatosis and Phospholipidosis Using Lipid Droplet-Selective Two-Photon Probes. <i>Analytical Chemistry</i> , 2020, 92, 11223-11231.	6.5	40
22	Design and synthesis of efficient heavy-atom-free photosensitizers for photodynamic therapy of cancer. <i>Chemical Communications</i> , 2020, 56, 11489-11492.	4.1	32
23	A fluorescent ESIPT-based benzimidazole platform for the ratiometric two-photon imaging of ONOO ⁻ <i>in vitro</i> and <i>ex vivo</i> . <i>Chemical Science</i> , 2020, 11, 7329-7334.	7.4	39
24	Two-photon imaging of hydrogen polysulfides in living cells and hippocampal tissues. <i>Sensors and Actuators B: Chemical</i> , 2020, 322, 128564.	7.8	29
25	Two-Photon and Multicolor Fluorogenic Bioorthogonal Probes Based on Tetrazine-Conjugated Naphthalene Fluorophores. <i>Bioconjugate Chemistry</i> , 2020, 31, 1545-1550.	3.6	15
26	Elevated TRPV4 Levels Contribute to Endothelial Damage and Scarring in Experimental Spinal Cord Injury. <i>Journal of Neuroscience</i> , 2020, 40, 1943-1955.	3.6	41
27	Discrimination between Human Colorectal Neoplasms with a Dual-Recognitive Two-Photon Probe. <i>Analytical Chemistry</i> , 2019, 91, 14705-14711.	6.5	19
28	Two-Photon Fluorescence Probe for Selective Monitoring of Superoxide in Live Cells and Tissues. <i>Analytical Chemistry</i> , 2019, 91, 14691-14696.	6.5	30
29	A two-photon fluorescent probe for colorimetric and ratiometric monitoring of mercury in live cells and tissues. <i>Chemical Communications</i> , 2019, 55, 1766-1769.	4.1	91
30	Ratiometric Detection of γ -Glutamyltransferase in Human Colon Cancer Tissues Using a Two-Photon Probe. <i>Analytical Chemistry</i> , 2019, 91, 9246-9250.	6.5	27
31	Unusual fluorescence of <i>o</i> -phenylazonaphthol derivatives with aggregation-induced emission and their use in two-photon cell imaging. <i>Chemical Communications</i> , 2019, 55, 6747-6750.	4.1	23
32	A Highly Sensitive Two-Photon Ratiometric Probe for Rapid Detection of the hNQO1 Enzyme in Colon Cancer Tissue. <i>Asian Journal of Organic Chemistry</i> , 2019, 8, 1707-1712.	2.7	9
33	Azulene-Derived Fluorescent Probe for Bioimaging: Detection of Reactive Oxygen and Nitrogen Species by Two-Photon Microscopy. <i>Journal of the American Chemical Society</i> , 2019, 141, 19389-19396.	13.7	125
34	A two-photon ratiometric probe for hydrogen polysulfide (H ₂ S _n): Increase in mitochondrial H ₂ S _n production in a Parkinson's disease model. <i>Sensors and Actuators B: Chemical</i> , 2019, 283, 810-819.	7.8	32
35	A two-photon ESIPT based fluorescence probe for specific detection of hypochlorite. <i>Dyes and Pigments</i> , 2018, 158, 526-532.	3.7	67
36	N-Heterocyclic Carbene Boranes as Reactive Oxygen Species-Responsive Materials: Application to the Two-Photon Imaging of Hypochlorous Acid in Living Cells and Tissues. <i>Angewandte Chemie - International Edition</i> , 2018, 57, 1567-1571.	13.8	127

#	ARTICLE	IF	CITATIONS
37	Pyrrrolidine dithiocarbamate reverses Bcl-xL-mediated apoptotic resistance to doxorubicin by inducing paraptosis. <i>Carcinogenesis</i> , 2018, 39, 458-470.	2.8	20
38	N-Heterocyclic Carbene Boranes as Reactive Oxygen Species-Responsive Materials: Application to the Two-Photon Imaging of Hypochlorous Acid in Living Cells and Tissues. <i>Angewandte Chemie</i> , 2018, 130, 1583-1587.	2.0	26
39	A Two-Photon Ratiometric Fluorescent Probe for Imaging of Hydrogen Peroxide Levels in Rat Organ Tissues. <i>ChemistryOpen</i> , 2018, 7, 53-56.	1.9	12
40	High-depth fluorescence imaging using a two-photon FRET system for mitochondrial pH in live cells and tissues. <i>Chemical Communications</i> , 2018, 54, 13531-13534.	4.1	48
41	Endoplasmic Reticulum-Targeted Ratiometric N-Heterocyclic Carbene Borane Probe for Two-Photon Microscopic Imaging of Hypochlorous Acid. <i>Analytical Chemistry</i> , 2018, 90, 12937-12943.	6.5	75
42	Near-IR Fluorescent Tracer for Glucose-Uptake Monitoring in Live Cells. <i>Bioconjugate Chemistry</i> , 2018, 29, 3394-3401.	3.6	22
43	A two-photon ratiometric probe for detection of hNQO1 enzyme activity in human colon tissue. <i>Sensors and Actuators B: Chemical</i> , 2018, 272, 203-210.	7.8	24
44	Carboxylate-Containing Two-Photon Probe for the Simultaneous Detection of Extra- and Intracellular pH Values in Colon Cancer Tissue. <i>Analytical Chemistry</i> , 2018, 90, 8058-8064.	6.5	18
45	A Two-Photon Fluorescent Probe for Imaging Endogenous ONOO ⁻ near NMDA Receptors in Neuronal Cells and Hippocampal Tissues. <i>Analytical Chemistry</i> , 2018, 90, 9347-9352.	6.5	71
46	Naphthalene-based fluorescent probes for glutathione and their applications in living cells and patients with sepsis. <i>Theranostics</i> , 2018, 8, 1411-1420.	10.0	31
47	Carboxylesterase-2-Selective Two-Photon Ratiometric Probe Reveals Decreased Carboxylesterase-2 Activity in Breast Cancer Cells. <i>Analytical Chemistry</i> , 2018, 90, 9465-9471.	6.5	49
48	Ratiometric Two-Photon Fluorescent Probe for Detecting and Imaging Hypochlorite. <i>Analytical Chemistry</i> , 2018, 90, 9510-9514.	6.5	86
49	Two-photon fluorescence sensors for imaging NMDA receptors and monitoring release of Zn ²⁺ from the presynaptic terminal. <i>Biosensors and Bioelectronics</i> , 2017, 91, 770-779.	10.1	24
50	Two-Photon Dye Cocktail for Dual-Color 3D Imaging of Pancreatic Beta and Alpha Cells in Live Islets. <i>Journal of the American Chemical Society</i> , 2017, 139, 3480-3487.	13.7	30
51	Visualization of vesicular transport from the endoplasmic reticulum to lysosome using an amidine derived two-photon probe. <i>Chemical Communications</i> , 2017, 53, 6097-6100.	4.1	14
52	Î±5-Syntrophin stabilizes catalase to reduce endogenous reactive oxygen species levels during myoblast differentiation. <i>FEBS Journal</i> , 2017, 284, 2052-2065.	4.7	8
53	An efficient two-photon fluorescent probe for human NAD(P)H:quinone oxidoreductase (hNQO1) detection and imaging in tumor cells. <i>Chemical Communications</i> , 2017, 53, 525-528.	4.1	56
54	A two-photon fluorescent probe for specific detection of hydrogen sulfide based on a familiar ESIPT fluorophore bearing AIE characteristics. <i>Chemical Communications</i> , 2017, 53, 4791-4794.	4.1	116

#	ARTICLE	IF	CITATIONS
55	Highly Selective and Sensitive Two-Photon Fluorescence Probe for Endogenous Peroxynitrite Detection and Its Applications in Living Cells and Tissues. <i>Analytical Chemistry</i> , 2017, 89, 8496-8500.	6.5	93
56	A ratiometric two-photon probe for Ca ²⁺ in live tissues and its application to spinal cord injury model. <i>Biomaterials</i> , 2017, 141, 251-259.	11.4	29
57	New Six-Membered pH-Insensitive Rhodamine Spirocyclic in Selective Sensing of Cu ²⁺ through C-C Bond Cleavage and Its Application in Cell Imaging. <i>ACS Omega</i> , 2017, 2, 8167-8176.	3.5	28
58	A viscosity sensitive fluorescent dye for real-time monitoring of mitochondria transport in neurons. <i>Biosensors and Bioelectronics</i> , 2016, 86, 885-891.	10.1	98
59	A Selective Imidazoline-2-thione-Bearing Two-Photon Fluorescent Probe for Hypochlorous Acid in Mitochondria. <i>Analytical Chemistry</i> , 2016, 88, 6615-6620.	6.5	160
60	A quadrupolar two-photon fluorescent probe for in vivo imaging of amyloid- β^2 plaques. <i>Chemical Science</i> , 2016, 7, 4600-4606.	7.4	49
61	A carboxylesterase-selective ratiometric fluorescent two-photon probe and its application to hepatocytes and liver tissues. <i>Chemical Science</i> , 2016, 7, 3703-3709.	7.4	54
62	Readily Accessible and Predictable Naphthalene-Based Two-Photon Fluorophore with Full Visible-Color Coverage. <i>Chemistry - A European Journal</i> , 2016, 22, 14166-14170.	3.3	10
63	Real-time monitoring of vesicle pH in an endocytic pathway using an EGF-conjugated two-photon probe. <i>Chemical Communications</i> , 2016, 52, 14007-14010.	4.1	14
64	Two-photon fluorescent probe for peroxynitrite. <i>Tetrahedron Letters</i> , 2016, 57, 715-718.	1.4	22
65	A ratiometric two-photon probe for quantitative imaging of mitochondrial pH values. <i>Chemical Science</i> , 2016, 7, 766-773.	7.4	118
66	One-Photon and Two-Photon Sensing of Biothiols Using a Bis-Pyrene-Cu(II) Ensemble and Its Application To Image GSH in the Cells and Tissues. <i>Analytical Chemistry</i> , 2015, 87, 3308-3313.	6.5	95
67	A cysteamine-selective two-photon fluorescent probe for ratiometric bioimaging. <i>Chemical Communications</i> , 2015, 51, 2407-2410.	4.1	34
68	Development of Imidazoline-2-thiones Based Two-Photon Fluorescence Probes for Imaging Hypochlorite Generation in a Co-Culture System. <i>Angewandte Chemie - International Edition</i> , 2015, 54, 4890-4894.	13.8	217
69	A Golgi-localized two-photon probe for imaging zinc ions. <i>Chemical Communications</i> , 2015, 51, 12099-12102.	4.1	42
70	Small-Molecule Two-Photon Probes for Bioimaging Applications. <i>Chemical Reviews</i> , 2015, 115, 5014-5055.	47.7	889
71	A hexaphenylbenzene based AIEE active two photon probe for the detection of hydrogen sulfide with tunable self-assembly in aqueous media and application in live cell imaging. <i>Chemical Communications</i> , 2015, 51, 15570-15573.	4.1	35
72	Mechanism of Cisplatin-Induced Cytotoxicity Is Correlated to Impaired Metabolism Due to Mitochondrial ROS Generation. <i>PLoS ONE</i> , 2015, 10, e0135083.	2.5	210

#	ARTICLE	IF	CITATIONS
73	Quinoline-Based Two-Photon Fluorescent Probe for Nitric Oxide in Live Cells and Tissues. <i>Analytical Chemistry</i> , 2014, 86, 308-311.	6.5	90
74	A small molecule two-photon fluorescent probe for intracellular sodium ions. <i>Chemical Communications</i> , 2014, 50, 1309-1312.	4.1	80
75	Ratiometric Two-Photon Fluorescent Probe for Quantitative Detection of β -Galactosidase Activity in Senescent Cells. <i>Analytical Chemistry</i> , 2014, 86, 10001-10005.	6.5	131
76	Red Emissive Two-Photon Probe for Real-Time Imaging of Mitochondria Trafficking. <i>Analytical Chemistry</i> , 2014, 86, 5638-5641.	6.5	62
77	Two-Photon Fluorescent Probes for Metal Ions in Live Tissues. <i>Inorganic Chemistry</i> , 2014, 53, 1794-1803.	4.0	72
78	Benzimidazole-Based Ratiometric Two-Photon Fluorescent Probes for Acidic pH in Live Cells and Tissues. <i>Journal of the American Chemical Society</i> , 2013, 135, 17969-17977.	13.7	306
79	A two-photon fluorescent probe for amyloid- β plaques in living mice. <i>Chemical Communications</i> , 2013, 49, 1303.	4.1	54
80	A Ratiometric Two-Photon Fluorescent Probe Reveals Reduction in Mitochondrial H_2O_2 Production in Parkinson's Disease Gene Knockout Astrocytes. <i>Journal of the American Chemical Society</i> , 2013, 135, 9915-9923.	13.7	383
81	Simultaneous Imaging of Mitochondria and Lysosomes by Using Two-Photon Fluorescent Probes for Zinc Ions and Thiols in Living Tissues. <i>Oxidative Medicine and Cellular Longevity</i> , 2013, 2013, 1-11.	4.0	20
82	Simultaneous Imaging of Mitochondria and Lysosomes by Using Two-Photon Fluorescent Probes. <i>Chemistry - A European Journal</i> , 2012, 18, 15246-15249.	3.3	42
83	A mitochondria-localized two-photon fluorescent probe for ratiometric imaging of hydrogen peroxide in live tissue. <i>Chemical Communications</i> , 2012, 48, 3518.	4.1	149
84	A Small Molecule Two-Photon Probe for Nitric Oxide in Living Tissues. <i>Chemistry - A European Journal</i> , 2012, 18, 12388-12394.	3.3	49
85	Dual-Color Imaging of Magnesium/Calcium Ion Activities with Two-Photon Fluorescent Probes. <i>Analytical Chemistry</i> , 2012, 84, 8110-8113.	6.5	52
86	Detection of Nickel in Fish Organs with a Two-Photon Fluorescent Probe. <i>Chemistry - A European Journal</i> , 2012, 18, 1953-1960.	3.3	26
87	Two-Photon Lyso trackers for in Vivo Imaging. <i>Journal of Organic Chemistry</i> , 2011, 76, 8113-8116.	3.2	44
88	A Mitochondrial-Targeted Two-Photon Probe for Zinc Ion. <i>Journal of the American Chemical Society</i> , 2011, 133, 5698-5700.	13.7	227
89	Ratiometric Detection of Mitochondrial Thiols with a Two-Photon Fluorescent Probe. <i>Journal of the American Chemical Society</i> , 2011, 133, 11132-11135.	13.7	348
90	Two-Photon Fluorescent Probes for Metal Ions. <i>Chemistry - an Asian Journal</i> , 2011, 6, 58-69.	3.3	127

#	ARTICLE	IF	CITATIONS
91	A Two-Photon Turn-On Probe for Lipid Rafts with Minimum Internalization. <i>ChemBioChem</i> , 2011, 12, 392-395.	2.6	22
92	Sodium-Ion-Selective Two-Photon Fluorescent Probe for In Vivo Imaging. <i>Angewandte Chemie - International Edition</i> , 2010, 49, 364-367.	13.8	101
93	Dual-Color Imaging of Sodium/Calcium Ion Activities with Two-Photon Fluorescent Probes. <i>Angewandte Chemie - International Edition</i> , 2010, 49, 6786-6789.	13.8	64
94	A Two-Photon Tracer for Glucose Uptake. <i>Angewandte Chemie - International Edition</i> , 2009, 48, 8027-8031.	13.8	55
95	Two-photon materials with large two-photon cross sections. Structure-property relationship. <i>Chemical Communications</i> , 2009, , 153-164.	4.1	102
96	Second-order nonlinear optical properties of octupolar molecules structure-property relationship. <i>Journal of Materials Chemistry</i> , 2009, 19, 7402.	6.7	69
97	Two-Photon Probes for Intracellular Free Metal Ions, Acidic Vesicles, And Lipid Rafts in Live Tissues. <i>Accounts of Chemical Research</i> , 2009, 42, 863-872.	15.6	530
98	Two-Photon Fluorescent Probes for Biomembrane Imaging: Effect of Chain Length. <i>ChemBioChem</i> , 2008, 9, 2830-2838.	2.6	32
99	Two-Photon Fluorescent Probes for Long-Term Imaging of Calcium Waves in Live Tissue. <i>Chemistry - A European Journal</i> , 2008, 14, 2075-2083.	3.3	50
100	Two-Photon Fluorescent Probes for Acidic Vesicles in Live Cells and Tissue. <i>Angewandte Chemie - International Edition</i> , 2008, 47, 2231-2234.	13.8	122
101	Two-Photon Fluorescent Probes for Intracellular Free Zinc Ions in Living Tissue. <i>Angewandte Chemie - International Edition</i> , 2008, 47, 5167-5170.	13.8	125
102	Two-Photon Absorption Properties of Alkynyl-Conjugated Pyrene Derivatives. <i>Journal of Organic Chemistry</i> , 2008, 73, 5127-5130.	3.2	102
103	Two-Photon Fluorescent Turn-On Probe for Lipid Rafts in Live Cell and Tissue. <i>Journal of the American Chemical Society</i> , 2008, 130, 4246-4247.	13.7	123
104	Magnesium Ion Selective Two-Photon Fluorescent Probe Based on a Benzo[h]chromene Derivative for in Vivo Imaging. <i>Journal of Organic Chemistry</i> , 2007, 72, 2088-2096.	3.2	136
105	A Two-Photon Fluorescent Probe for Lipid Raft Imaging: C-Laurdan. <i>ChemBioChem</i> , 2007, 8, 553-559.	2.6	228
106	Environment-Sensitive Two-Photon Probe for Intracellular Free Magnesium Ions in Live Tissue. <i>Angewandte Chemie - International Edition</i> , 2007, 46, 3460-3463.	13.8	151
107	A Two-Photon Fluorescent Probe for Calcium Waves in Living Tissue. <i>Angewandte Chemie - International Edition</i> , 2007, 46, 7445-7448.	13.8	102
108	Design of molecular two-photon probes for in vivo imaging. 2H-Benzo[h]chromene-2-one derivatives. <i>Tetrahedron Letters</i> , 2007, 48, 2791-2795.	1.4	56

#	ARTICLE	IF	CITATIONS
109	Ketene-forming elimination reactions from aryl phenylacetates promoted by R_2NH in MeCN: effects of base-solvent and σ -phenyl group. <i>Journal of Physical Organic Chemistry</i> , 2007, 20, 685-689.	1.9	5
110	First hyperpolarizabilities of hexa(ethynyl)benzene derivatives: effect of conjugation length. <i>Journal of Materials Chemistry</i> , 2006, 16, 2273.	6.7	24
111	Metal Ion Sensing Novel Calix[4]crown Fluoroionophore with a Two-Photon Absorption Property. <i>Journal of Organic Chemistry</i> , 2006, 71, 8016-8022.	3.2	71
112	First Hyperpolarizabilities of 1,3,5-Tricyanobenzene Derivatives: Origin of Larger β^2 Values for the Octupoles than for the Dipoles. <i>ChemPhysChem</i> , 2006, 7, 206-212.	2.1	77
113	Molecular two-photon sensor for metal ions derived from bis(2-pyridyl)amine. <i>Chemical Physics Letters</i> , 2005, 410, 312-315.	2.6	35
114	Two-Photon Absorption Properties of 2,6-Bis(styryl)anthracene Derivatives: Effects of Donor-Acceptor Substituents and the π Center. <i>Chemistry - A European Journal</i> , 2005, 11, 4191-4198.	3.3	75
115	Two-Photon Dyes Containing Heterocyclic Rings with Enhanced Photostability. <i>Chemistry - A European Journal</i> , 2005, 11, 6386-6391.	3.3	32
116	Two-Photon Sensor for Metal Ions Derived from Azacrown Ether. <i>Journal of Organic Chemistry</i> , 2004, 69, 5749-5751.	3.2	73